


Varomedshoes: which size fits?

As a matter of principle: all article numbers, starting with the same double-digit number (such as 60.xxx-xx), have the same length and width.



How to do it:

1. measure the footlength at its longest part
2. add 5 - 7 mm, the toe should have enough space
3. look for the next bigger size in the chart of your article group, column a

Article group 58

sizes:	36	37	38	39	40	41	42	43	44	45	46	47/48
a	231	238	245	251	258	265	271	278	285	291	298	310
b	83	85	87	89	91	93	95	97	99	101	103	107
c	60	61,5	63	64,5	66	67,5	69	70,5	72	73,5	75	77

All measurements are shown in millimeters.

Article groups 60, 75, 76, 77, 78

Größen:	36	37	38	39	40	41	42	43	44	45	46	47/48
a	239	246	253	259	266	273	279	286	293	299	306	313
b	91	93	95	97	99	101	103	105	107	109	111	113
c	60,5	62	63,5	65	66,5	68	69,5	71	72,5	74	75,5	77

All measurements are shown in millimeters.